
Perasafe 교육자료 II

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"Antibiotics" are strictly pharmaceuticals of natural origin which kill or inhibit bacteria. They have a very specific mode of action and are of low enough toxicity to be swallowed or injected into the body. The term "antibacterial" was originally used for similar pharmaceuticals made by chemical synthesis. Today, however, the term is applied to hygiene products which kill or inhibit bacteria outside the body, but the ingredients are very different and have a general rather than specific mode of action.

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chlorinated derivatives of phenols이며, 이것은

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


Benzalkonium chloride


Cetrimide

Quaternary ammonium compounds with just one "fatty" chain, like the two above, are strongly biocidal. As more fatty chains are added, the molecule loses its biocidal activity, but gains surfactant activity. Such compounds are used in fabric conditioners and in specialised detergents.

Did you know that as well as killing germs to prevent disease, or preserving food, "biocides" that kill or prevent the growth of micro-organisms are used in all kinds of products - paints, plastics, wallpaper, cutting oils, wood, textiles and leather for example - to prevent spoilage and keep them safe to use? Biocides are also used to prevent microbes "fouling" pipes and equipment in factories that use a lot of water - paper mills and power stations for example.

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INFORMATION BOX - DISINFECTANTS AND THEIR ORIGINS

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Standards can always be improved

| | 2% GLUTARALDEHYDE | 1.62% PERASAFE |
|----------------------------|------------------------------|---------------------------|
| STERILISATION | 10 hours | 10 minutes |
| DISINFECTION | 10 minutes | 5 minutes |
| MYCOBACTERIA | 60 minutes | 10 minutes |
| SPORES | 30 minutes | 1 minute |
| PS AERUGINOSA | 10 minutes | 5 minutes |
| ENVIRONMENT | Damaging | Friendly |
| VENTILATION COST | High | None |
| TOXIC FUMES | High | None |
| PROTECTIVE CLOTHING | Required | None |
| REQUIRES ACTIVATION | Yes | Yes |
| ENZYMATIC PRE-CLEAN | Recommended | Not required |

Disinfection vs Sterilisation

□ Disinfection

“a process giving a **5-log reduction (99.999%)** kill of viruses, bacteria and fungi”

□ Sterilisation

“a process removing all trace of living material”

- Sterilisation tests call for spore kills of the order of **12-log reduction, i.e. 99.999999999999% kill**

Disinfection vs Sterilisation

❖ When to disinfect vs When to sterilise

❑ IDEALLY, all instruments and equipment should be ‘cleaned’ and then sterilised, preferably by autoclave (steam + pressure)

❑ However

- Autoclaves are simply not available in all markets
- Those which are available may not be well-maintained or reliable
- Key instruments (flexible endoscopes) are heat labile and cannot be autoclaved
- Sterilization process techniques (autoclaving, irradiation, ETO, etc) are time-consuming – as is aldehyde sterilization – and may be prohibitively expensive

Disinfection vs Sterilisation

❖ Key reasons why PeraSafe is needed

- ❑ All Instruments should be disinfected and sterilised
- ❑ Heat-labile endoscopes require rapid sterilization/disinfection WITHOUT the toxic effects of the aldehydes
- ❑ A safe, fast acting replacement for the aldehydes is required in those countries dependent on “cold liquid sterilants”

Disinfection vs Sterilization

❖ The Medical Requirements of an Instrument disinfectant/sterilant

The essential properties –

- ❑ Total biocidal spectrum
 - Producing complete inactivation of spores mycobacteria and viruses on short contact times

- ❑ Rapid in action
 - Permitting rapid disinfection and sterilisation between patients which is essential due to the limited number of instruments available

Disinfection vs Sterilisation


❖ The Medical Requirements of an Instrument disinfectant/sterilant

The essential properties –

- ❑ Total biocidal spectrum
 - Producing complete inactivation of spores mycobacteria and viruses on short contact times

- ❑ Rapid in action
 - Permitting rapid disinfection and sterilisation between patients which is essential due to the limited number of instruments available

Simple to activate,
forms peracetyl ions at pH 8

The Ultimate General Purpose Disinfectant

The 10 Minute Medical Instrument High Level Disinfectant and Sterilant